



OMICS and Precision Medicine - Full Workflows for Lipidomics and Metabolic Profiling

Jeremiah D. Tipton, Ph.D.
Account Manager TN & MS

The Promise of Quantitative Systems Biology

- **Human Health**

- Biomarkers will help monitor wellness, detect early disease, monitor disease progression, and stratify patients for more tailored care.



- **Agriculture and Food**

- Bioengineering of key crops will enable the optimization of production, nutrition, taste, durability, and resistance to infectious agents.





- **Energy and Environment**

- Bioengineering based on a systems level knowledge for the development of alternative bioenergy sources and improved carbon sequestration.




The Challenges in Systems Biology




3 © 2016 AB Sciex



The Complexity of Systems Biology for Personalized | Precision Medicine




OneOmics™ Project



OMICS Technologies


The OMICS Spectrum



Wenk et al. Nature 2005

OMICS Technologies

The OMICS Spectrum










Wenk et al. Nature 2005

7 © 2016 AB Sciex



SCIEX Metabolomics Portfolio

Discovery	Translation/Validation	Clinical Utilization
 <p>Global Untargeted Screening</p> <ul style="list-style-type: none"> • 1000's of <i>features</i> in 10's of samples • SWATH® / IDA Acquisition • Brings data analysis from weeks to days  <p>Global Targeted Screening</p> <ul style="list-style-type: none"> • Profile over 500 known <i>metabolites</i> in 10's of samples • SWATH® / IDA Acquisition • SCIEX Accurate Mass Spectral Library: processing to visualization in minutes 	 <p>Targeted Profiling</p> <ul style="list-style-type: none"> • 100s of putative biomarkers on 100-1000s of samples • QTRAP® Systems: Industry standard robustness • Kits / reagents / existing methods: <ul style="list-style-type: none"> • aTRAQ™ Reagent for AAA • Ampliflex reagents for enhanced sensitivity • MRM assays for Vitamin D, Hormones/Steroids 	 <p>Validated Assays</p> <ul style="list-style-type: none"> • 10s of biomarkers and 100,000s of samples/yr. • SCIEX Triple Quad™ 4500 MD • MultiQuant™ Software MD • Cliquant™ Software MD 

Putative biomarkers

Samples

Applications FluxOMICS and Pathways

Targeted Metabolomics - FluxOMICS

Quantitative and Qualitative Metabolomics for the Investigation of Intracellular Metabolism

Targeted Analysis on the QTRAP® 5500 System and Reverse-Phase Ion-Pairing Chromatography

Douglas McCloskey¹ and Baljit K. Ubhi²

¹Department of Bioengineering, University of California, San Diego, CA, USA, ²SCIEX, USA

Liquid chromatography-mass spectrometry (LC-MS) provides a powerful analytical tool for understanding and monitoring intracellular metabolism by measuring the metabolome. The study of intracellular metabolism of model organisms, such as *E. coli*, is vital to further our biochemical knowledge¹, to develop new pharmaceuticals that target harmful pathogens², and to improve industrial applications that aim to metabolically engineer bacteria in order to produce commodity chemicals from renewable resources³. Paramount to these endeavors is the ability to reliably and accurately measure the intracellular metabolome. For microorganisms, the compounds of most interest comprise intermediates of high flux pathways such as glycolysis, the pentose phosphate pathway, the citric acid cycle, amino acid metabolism, as well as energy and redox cofactors such as ATP and NADH (Figure 1). By measuring the absolute metabolite levels of such compounds, one is able to calculate reaction and pathway thermodynamics⁴ and infer *in vivo* enzyme kinetics⁵. In addition, when microorganisms are grown on a specifically chosen labeled substrate (e.g. 1-¹³C glucose) during a metabolic labeling experiment, the isotopomer distribution of intracellular compounds can be used to calculate the absolute flux through specific reactions of interest⁶.



In this work, the QTRAP® 5500 system (a hybrid triple quadrupole linear ion trap mass spectrometer) was used to implement both quantitative and qualitative workflows aimed at measuring anionic and polar compounds of intracellular metabolism.

Key Features of the QTRAP® 5500 System for Qualitative and Quantitative Metabolomics

Targeted Metabolomics - FluxOMICS




Figure 1. Metabolomics of Intracellular Pathways for Investigations into the Biochemistry of Microorganisms. High flux pathways such as those shown are key to generating a biological picture and targeted metabolomics strategies provide a robust quantitative strategy for monitoring changes.

11 © 2016 AB Sciex






Targeted Metabolomics - FluxOMICS

Metabolite	Q1	Q3	RT	DP	EP	CE	CXP	MS Window	Quantifier
L-Arginine-UC13	179	136	0.7	-125.3	-10	-19.1	-6.9	30	1
L-Arginine	173	131	0.7	-125.3	-10	-19.1	-6.9	30	1
L-Ornithine-UC13	136	89	0.7	-61.8	-10	-17.4	-7.9	30	1
L-Ornithine	131	85	0.7	-61.8	-10	-17.4	-7.9	30	1
L-Asparagine-UC13	135	118	0.92	-89.7	-10	-16.1	-7	30	1
L-Asparagine	131	114	0.92	-89.7	-10	-16.1	-7	30	1
L-Asparagine	131	95	0.92	-51	-10	-17.1	-11.3	30	2
L-Serine-UC13	107	76	0.92	-72.8	-10	-15.9	-8.7	30	1
L-Serine	104	74	0.92	-72.8	-10	-15.9	-8.7	30	1
L-Serine	104	42	0.92	-73.9	-10	-18.7	-12.2	30	2
L-Citrulline-UC13	180	136	0.93	-76.4	-10	-19.1	-6.4	30	1
L-Citrulline	174	131	0.93	-76.4	-10	-19.1	-6.4	30	1
L-Glutamine-UC13	150	132	0.93	-19.9	-10	-13.3	-17.7	30	1
L-Glutamine	145	127	0.93	-19.9	-10	-13.3	-17.7	30	1
L-Glutamine	145	101	0.93	-20.2	-10	-14.3	-6.8	30	2
Hexose_Pool_fru_glc-D-UC13	185	61	0.94	-36.2	-10	-23.1	-5.8	30	1
Hexose_Pool_fru_glc-D	179	59	0.94	-36.2	-10	-23.1	-5.8	30	1
Hexose_Pool_fru_glc-D	179	89	0.94	-37.4	-10	-12.9	-6.7	30	2
L-Cystine-UC13	245	123	0.94	-59.7	-10	-17	-14.8	30	1
L-Cystine	239	120	0.94	-59.7	-10	-17	-14.8	30	1
L-Threonine-UC13	122	76	0.94	-46.5	-10	-14.2	-6.7	30	1
L-Threonine	118	74	0.94	-46.5	-10	-14.2	-6.7	30	1
L-Alanine-UC13	91	46	0.95	-14.2	-10	-19.9	-14	30	1
L-Alanine	88	45	0.95	-14.2	-10	-19.9	-14	30	1
L-Alanine	88	42	0.95	-11.2	-10	-19.2	-17	30	2
Cytidine-UC13	251	113	1.15	-100	-4	-16.1	-5.2	30	1
Cytidine	242	109	1.15	-100	-4	-16.1	-5.2	30	1
L-Histidine-UC13	160	143	1.15	-56.7	-10	-18.8	-14.8	30	1
L-Histidine	154	137	1.15	-56.7	-10	-18.8	-14.8	30	1
L-Histidine	154	93	1.15	-57.6	-10	-22.8	-9	30	2
Uracil-UC13	115	43	1.15	-100	-10	-26	-7	30	1
Uracil	111	42	1.15	-100	-10	-26	-7	30	1
5-Aminoimidazole-4-carboxamide riboside-UC13	266	129	1.2	-100	-10	-22	-9	30	1

12 © 2016 AB Sciex



Targeted Metabolomics - FluxOMICS



13 © 2016 AB Sciex

Targeted Metabolomics - FluxOMICS




Figure 3. Simultaneous Quantitative Analysis with Qualitative Confirmation of L-Glutamate using the QTRAP® 5500 System. A) The primary and secondary transitions (blue and green) for glu-L are monitored, along with the uniformly heavy carbon labeled analog (red). A full scan MS/MS spectrum (EPI scan) is triggered when the primary and secondary transitions reach a predefined threshold. B) To confirm the identity of the MRM signal, the MS/MS was matched with a greater than 90% match to the reference spectrum (taken for pure standards for glutamate).

14 © 2016 AB Sciex




Figure 4. Qualitative Method for Characterizing the Isotopomers. Enabled by the QTRAP 5500® System. The advantage of the qualitative method is highlighted here in this example measuring unlabeled and fully labeled malic acid (mal-L) in *E. coli*. Several transitions corresponding to the isotopomer distribution are monitored per compound. An Enhanced Resolution (ER) scan is triggered when one of the isotopomer transitions reaches a predefined threshold and provides MS and isotope ratio information. An EPI scan is also triggered and provides information regarding the location of the heavy label (if present).

Targeted Metabolomics - FluxOMICS




Figure 6. Heat Map of Absolute Metabolite Level Differences Between *E. coli* Grown in Minimal Media and Complex Media. The top 25 most significant metabolites (x-axis) are clustered (Pearson's R, complete linkage) according to their significance and fold-change (Student's t-test). The media conditions (y-axis) are clustered (Pearson's R, complete linkage) according to similarity in metabolite level of the top 25 most significant metabolite changes.

15 © 2016 AB Sciex



The Lipidizer™ Platform: A Revolutionary Tool for Understanding the Role of Lipids in Disease

For Research Use Only. Not for use in diagnostic procedures. RUO-MKT-11-4185-A


Industrializing Lipidomics

Lipidyzer™ Platform

- Robust MS System
- Class separation using DMS
- Unique internal standard strategy with validated kits
- Data Visualization including heat maps, box & whisker plots, etc



Lipidyzer™ Platform



17 © 2016 AB Sciex



Demonstrating the Power of Lipidyzer™ Platform

Benefits

- Specificity – Differential Mobility Spectrometry (DMS)
 - Eliminating Quantitative Bias – Novel Internal Standards
 - Coverage – Complex Lipid Metabolism
 - Sensitivity and Precision - Assay
 - Robustness – Assay and Platform
 - Ease of Use - Platform
-
- Samples Sets used for Biological Validation

18 © 2016 AB Sciex



Complex Lipids are like a Matrix

- Lipid are present in classes that have concentrations and compositions (important for level of metabolism)
 - Concentration = sum of the FAs for any given class (column)
 - Composition = relative abundances of each FA (or species) across many classes (rows)

FATTY ACIDS	LIPID CLASSES							
	CE	TAG	DAG	FFA	PC	PE	LPC	LPE
14:0								
16:0								
18:0								
20:0								
24:0								
14:1								
16:1								
18:1								
20:1								
18:2								
18:3								
20:2								
20:3								
20:4								
20:5								
22:4								
22:5								
22:6								

Sum = concentration

Sum = composition

19 © 2016 AB Sciex



Complex Lipids are like a Matrix

- Lipid are present in classes that have concentrations and compositions (important for level of metabolism)
 - Concentration = sum of the FAs for any given class (column)
 - Composition = relative abundances of each FA (or species) across many classes (rows)
- When FA metabolism is altered there is the ability to change FA composition of all classes

FATTY ACIDS	LIPID CLASSES							
	CE	TAG	DAG	FFA	PC	PE	LPC	LPE
14:0								
16:0								
18:0								
20:0								
24:0								
14:1								
16:1								
18:1								
20:1								
18:2								
18:3								
20:2								
20:3								
20:4								
20:5								
22:4								
22:5								
22:6								

20 © 2016 AB Sciex



Complex Lipids are like a Matrix

- Lipid are present in classes that have concentrations and compositions (important for level of metabolism)
 - Concentration = sum of the FAs for any given class (column)
 - Composition = relative abundances of each FA (or species) across many classes (rows)
- When FA metabolism is altered there is the ability to change FA composition of all classes
- When lipid class metabolism is altered there is the ability to change all members of the class

		LIPID CLASSES								
		CE	TAG	DAG	FFA	PC	PE	LPC	LPE	
FATTY ACIDS	14:0									
	16:0									
	18:0									
	20:0									
	24:0									
	14:1									
	16:1									
	18:1									
	20:1									
	18:2									
	18:3									
	20:2									
	20:3									
	20:4									
	20:5									
	22:4									
	22:5									
	22:6									



21 © 2016 AB Sciex

What is needed from a Lipid Platform

		LIPID CLASSES								
		CE	TAG	DAG	FFA	PC	PE	LPC	LPE	
FATTY ACIDS	14:0									
	16:0									
	18:0									
	20:0									
	24:0									
	14:1									
	16:1									
	18:1									
	20:1									
	18:2									
	18:3									
	20:2									
	20:3									
	20:4									
	20:5									
	22:4									
	22:5									
	22:6									

1) Specificity

- A non-specific method (e.g. PC 36:2) does not allow mapping to the elements of the matrix

2) Quantitation

- A non-quantitative approach does not allow accurate summing of the rows and columns

3) Comprehensive Coverage













- A partially complete matrix is difficult to interpret

22 © 2016 AB Sciex



Broad Range of IS to Normalize Quantitative Data

- Diversity of fatty acid chain lengths and degrees of unsaturation result in differential fragmentation efficiency which impacts quantitation
- Multiple IS that reflect the diversity of lipid molecular species

PHOSPHATIDYLCHOLINE (PC) INTERNAL STANDARD MIX				
	STRUCTURE	FATTY ACID	POS	%
		FA16:1 - Palmitoleic acid	sn-2	5
		FA18:1 - Oleic acid	sn-2	20
		FA18:2 - Linoleic acid	sn-2	20
		FA18:3 - α -Linolenic acid	sn-2	5
		FA20:3 - Dihomo- γ -linolenic acid	sn-2	5
		FA20:4 - Arachidonic acid	sn-2	20
		FA20:5 - Eicosapentaenoic acid	sn-2	5
		FA22:4 - Eicosatetraenoic acid	sn-2	5
		FA22:5 - Docosapentaenoic acid	sn-2	5
		FA22:6 - Docosahexaenoic acid	sn-2	10
	d916:0 - Labeled palmitic acid	sn-1	100	

- Each lipid class has multiple internal standards at concentrations that reflect those found in biology

23 © 2016 AB Sciex




Specificity Offered by SelexION® Technology

Removal of Isobaric Interferences

- Problem: The Q1 isolation window during MS/MS is ~ 1.2 Da, which increases number of potential isobars

Ionogram: Separation by Lipid Head Group in the Gas Phase




Differential Mobility Spectrometry-Driven Shotgun Lipidomics *Anal. Chem.* 2014, 86, 9662-9669 10.1021/ac5021744

24 © 2016 AB Sciex



How Does SelexION™ Technology Separate Ions?

Differential Mobility Spectrometry (DMS) separates molecules using planar geometry




25 © 2016 AB Sciex



Separation of Lipid Classes Using SelexION™ Technology

Effects of Modifier on COV Values




26 © 2016 AB Sciex



The Lipidzyzer Eliminates Quantitative Bias

Multiple internal standards per class provide accurate quantitation



27 © 2016 AB Sciex



Full Coverage of Complex Lipid Metabolism

Coverage and Depth

- Over 1100 molecular species across 13 lipid classes
- Lipidzyzer™ Platform provides 6 measurements:
 1. Lipid Class Concentration
 2. Lipid Species Concentration
 3. Fatty Acid Concentration
 4. Lipid Class Composition
 5. Lipid Species Composition
 6. Fatty Acid Composition

Fraction	Lipid Classes	Number of Species*
Neutral Lipids	Triacylglycerols (TAG)	502
	Diacylglycerols (DAG)	67
	Free Fatty Acids (FFA)	28
	Cholesterol Esters (CE)	34
Polar Lipids	Phosphatidylcholines (PC)	161
	Phosphatidylethanolamines (PE)	233
	Lysophosphatidylcholines (LPC)	28
	Lysophosphatidylethanolamines (LPE)	28
	Sphingomyelins (SM)	16
	Ceramides (CER)	56

*The Ceramides listed above includes the further three classes, DCER, HCER and LCER


28 © 2016 AB Sciex



Lipidizer™ Platform Sensitivity and Precision

5 Day Study

- Validated across 4 instruments and 3 labs the instruments detected similar numbers of lipid species and with similar precision
- >675 lipid species with RSD <20% in this control sample




29 © 2016 AB Sciex



Lipidizer™ Platform Robustness

Precision: 6 Day Study in Plasma (Total Class Concentration)

- 11/13 Classes <10% RSD over 6 days



Lipid Class	%RSD
Cholesterol esters (CE)	4.7
Ceramides (CER)	6.4
Diacylglycerols (DAG)	9.1
Free Fatty Acids (FFA)	4.5
Hexosylceramides (HCER)	5.8
Lactosylceramides (LCER)	6.9
Lysophosphocholines (LPC)	9.9
Lysophosphoethanolamines (LPE)	10.9
Phosphocholines (PC)	5.1
Phosphoethanolamines (PE)	8.9
Sphingomyelins (SM)	5.4
Triacylglycerols (TAG)	5.1
Dihydroceramides (DCER)	29.5*

*Note: DCERs are present at exceedingly low levels in plasma

30 © 2016 AB Sciex




Quickest Route to Success

Significant changes between Young, Old and Calorie Restricted

- Univariate statistical analysis
 - Non-parametric Wilcoxon t test – FDR corrected q value < 0.05
 - Fold change > 1.5

Lipidyzer™ demonstrated:

- Fastest route to successful data ready for interpretation
- Larger number of detected species
- Quantitative data on all species detected
- Allows mapping data to biochemical pathways



37 © 2016 AB Sciex

Conclusions

Demonstrating the Power of Lipidyzer™ Platform

- Specificity – Differential Mobility Spectrometry (DMS)
- Eliminating Quantitative Bias – Novel Internal Standards
- Coverage – Over 1100 Molecular Species across 13 lipid classes
- Sensitivity and Precision – Quantitate ~700 species <20%CVs
- Robustness – <10%CVs over 6 days
- Ease of Use – Lipidomics Workflow Manager

- Biological Validation of the Lipidyzer Platform
 - Preeclampsia Pilot Study
 - Hypertriglyceridemia Study
 - Weight Loss Study
 - Aging Pilot Study

38 © 2016 AB Sciex



Acknowledgements

SCIEX

- Pauline Vollmerhaus
- Eva Duchoslav
- Paul Baker
- Scott Hamilton
- Marjorie Minkoff
- Larry Campbell
- Yves Le Blanc
- Aaron Hudson
- Mark Cafazzo

METABOLON

- Steve Watkins
- Alex Conner
- Annie Evans
- Richard Robinson
- UyenThao Nguyen
- Joseph Yu
- Luke Miller
- Sarada Tanikella
- Peter Zhu
- Corey DeHaven

Avanti Polar Lipids

- Walt Shaw
- Lisa Connelly

Collaborators

- Katherine Williams & Susan Fisher (UCSF)
- Andy Dannenberg (Weill College Medicine)
- Kevin Contrepois & Mike Snyder (Stanford)

Legal Acknowledgements:

For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license. RUO-MKT-11-4185-A

© 2016 AB SCIEX.

39 © 2016 AB Sciex

