

UAB CT Scanners

Before placing a patient on the scanner, ensure you have an accurate weight. No patient exceeding the below limits should be scanned. To do so would have a high risk of causing serious scanner damage, leading to possible denial of or delay in service to patients.

Vendor	Model	Table Weight
The Kirklint Clinic (TKC)		
Siemens	Force	500 lbs
GE	Discovery	400 lbs
Highlands		
Philips	i256	650 lbs
Philips	Brilliance 16	300 lbs
Gardendale		
GE	Revolution 128	500 lbs
Leeds		
GE	Revolution 256	500 lbs
Acton Road		
Philips	16 (iPatient)	650lbs
UAB Hospital (NP)		
Siemens	Edge	500lbs
Siemens	Edge Bariatric top	676 lbs
Siemens	Force	500 lbs
Siemens	Force Bariatric top	676 lbs
Siemens	Alpha	676 lbs
Philips	i256	650 lbs