

UAB MRI Scanners

Before placing a patient on the scanner, ensure you have an accurate weight. No patient exceeding the below limits should be scanned. To do so would have a high risk of causing serious scanner damage, leading to possible denial of or delay in service to patients.

Vendor	Model	Table Weight
The Kirklint Clinic (TKC)		
Philips	Achieva	550 lbs
Philips	Achieva	550 lbs
Siemens	Skyra	550 lbs
GE	Twin Speed	350 lbs
Philips	Ingenia	550 lbs
Highlands		
Siemens	Symphony	550 lbs
Siemens	Prisma	550 lbs
Siemens	FreeMax	700 lbs
GE	Premier	550 lbs
Gardendale		
GE	Artist	500 lbs
Leeds		
GE	450W	500 lbs
Advanced Imaging Facility (AIF)		
GE	Signa PET-MR	500 lbs
UAB Hospital (NP)		
Philips	Ingenia	550 lbs
Philips	Ingenia	550 lbs
Siemens	Sola	
Siemens	Area	550lbs