

## UAB CT Scanners

Before placing a patient on the scanner, ensure you have an accurate weight. No patient exceeding the below limits should be scanned. To do so would have a high risk of causing serious scanner damage, leading to possible denial of or delay in service to patients.

Vendor	Model	Table Weight
<b>The Kirklin Clinic (TKC)</b>		
Siemens	Force	500 lbs
GE	Discovery	400 lbs
<b>Highlands</b>		
Philips	i256	650 lbs
Canon	Aquilion ONE/INSIGHT	660 lbs
<b>Gardendale</b>		
GE	Revolution 128	500 lbs
<b>Leeds</b>		
GE	Revolution 256	500 lbs
<b>Acton Road</b>		
Philips	iPatient 16	650 lbs
<b>UAB Hospital (NP)</b>		
Siemens	Edge	500 lbs
Siemens	Edge Bariatric Top	676 lbs
Siemens	Force	500 lbs
Siemens	Force Bariatric Top	676 lbs
Siemens	Alpha	676 lbs
Philips	i256	650 lbs